

Bed Pan Rack with Dip Tray



- Stainless Steel Construction
- Bottle rack
- Stainless Steel Backing Plate
- Pan rack
- Removable Dip Tray
- Total depth 270mm

Warranty Information

- 5 Years

 Contract approved

Support @ Home

Function: Collecting urine and faeces

Category: Prescribed

Click [here](#) for more information

Compliance

ISO 15883/AS 4187:2019

Code	Display Name	Total Length	Total Width	Total Height
UPA200276	4 Pan / 4 Bottle - Stainless Steel	270 mm	520 mm	650 mm
UPA200277	6 Pan / 6 Bottle - Stainless Steel	270 mm	800 mm	650 mm
UPA200278	8 Pan / 8 Bottle - Stainless Steel	270 mm	1070 mm	650 mm