

Weigh Chair Digital



- Ideal for accurate and convenient patient weighing
- Easy to read weight display
- Swivel wheels with safety brakes plus built in footrests
- Easy to clean surfaces essential in health care environments
- 100 hours use from regular batteries



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load
LSS377700	Weigh Chair Digital	3469	1	635 mm	680 mm	880 mm	150 kg