

K Care Hand Made Commode Seat - Closed



- Superior quality white marine grade vinyl provides durability and comfort
- Vinyl is wrapped, stapled and sealed, making cleaning easier with no seam on the inside of the aperture
- Comes with plastic seat runners for easy attachment and removal for cleaning
- Seat height is 55 mm



Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year



Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width
BTC084900	Closed Front Commode Seat - Hand Made - 320 mm Wide	430 mm	320 mm
BTC073800	Closed Front Commode Seat - Hand Made - 400 mm Wide	430 mm	400 mm
BTC073700	Closed Front Commode Seat - Hand Made - 440 mm Wide	430 mm	440 mm
BTC062400	Closed Front Commode Seat - Hand Made - 500 mm Wide	430 mm	500 mm
BTC062500	Closed Front Commode Seat - Hand Made - 550 mm Wide	430 mm	550 mm
BTC062600	Closed Front Commode Seat - Hand Made - 600 mm Wide	430 mm	600 mm
BTC062700	Closed Front Commode Seat - Hand Made - 650 mm Wide	430 mm	650mm

Code	Display Name	Seat Depth	Seat Width
------	--------------	------------	------------