

Aspire Montana HD Lift Recline Chair



- True Bariatric Seating Solution.
- Dual Action: Unlike the standard Montana which is a single action, this chair is available in a dual action configuration with independent control of the back, leg rest, and lift movements
- Interchangeable Backrests: Three configurations that can easily be changed due to their velcro fastening with standard (pictured), waterfall, and reinforced lateral support
- Custom Sizing: Capability to offer clients a custom-sized lift chair with the ability to specify the seat width, depth, height and armrest height. This option is available via the custom order form
- Head Cushion: Included as standard (if chosen) for extra comfort and support
- Waterfall and Immersion Backrests are custom and need to be ordered
- Legrest capacity
 - Montana HD (CHP227405): 160 kg
 - Montana HD450 (CHP227410): 140 kg
- Looking to protect your Electric Lift Recline Chair? Try our [Universal Chair Protector](#)

 Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Warranty Information

- 2 Years Mechanism & Electronics
- 2 Years Upholstery
- Limited Lifetime Frame

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Backrest Height
CHP227405	Aspire Montana HD Lift Recline Chair - 254kg	900 mm	1270 mm	254 kg	500 mm	650 mm	810 mm
CHP227410	Aspire Montana HD450 Lift Recline Chair - 450kg	1100 mm	1270 mm	450 kg	500 mm	850 mm	810 mm