

## Jobst Elvarex Handpiece w/ Full Fingers



### Jobst Elvarex Armwear Overview:

- Jobst Elvarex garments feature a special low-stretch, double-knit fabric, which exerts a light massaging effect on the initial lymphatics, and minimizes re-oedematisation
- The high cotton content, flat knit construction and short-stretch thread means that Elvarex products provide an advanced anatomical fit with maximum comfort
- Compression: Class 1: 20 - 30mmHg

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information



Code	Display Name	Calf/Arm Circumference	Elbow Circumference	Wrist Circumference
PTS642173	JOBST Elvarex - Handpiece with full fingers (Right) - 20-30mmHg - Sml	230 - 260 mm	240 - 260 mm	160 - 180 mm
PTS642174	JOBST Elvarex - Handpiece with full fingers (Left) - 20-30mmHg - Sml	230 - 260 mm	240 - 260 mm	160 - 180 mm
PTS642175	JOBST Elvarex - Handpiece with full fingers (Right) - 20-30mmHg - Med	260 - 280 mm	260 - 290 mm	180 - 200 mm
PTS642176	JOBST Elvarex - Handpiece with full fingers (Left) - 20-30mmHg - Med	260 - 280 mm	260 - 290 mm	180 - 200 mm
PTS642177	JOBST Elvarex - Handpiece with full fingers (Right) - 20-30mmHg - Lge	280 - 300 mm	290 - 320 mm	200 - 230 mm
PTS642178	JOBST Elvarex - Handpiece with full fingers (Left) - 20-30mmHg - Lge	280 - 300 mm	290 - 320 mm	200 - 230 mm