

Kerry Homecare Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- Versatile unit that dismantles for transport and storage, suitable for domestic use
- Epoxy powder coated frame
- Single speed handset control
- Emergency lowering
- Wide push handle for excellent manoeuvrability
- 125 mm rear total lock castors and 75 mm front twin swivel castors
- Manual foot operated leg gating for close access
- Suits a large range of patient slings
- Supplied with battery and charger



Warranty Information

- 2 Years