

## The Emergency Cart



- Adjustable push handle for quick and easy operation to enable users to set push handle at the correct ergonomic height
- Utility hook rail
- Adjustable intravenous pole
- Slide out shelf fitted with runners for smooth, stable operation
- Raised drawer handles designed for grip and easy cleaning
- 5 drawers fitted with soft close runners and damper - prevents drawers from opening while in motion
- Corner buffers
- Directional lock castor
- Diameter 125 anti-static total lock castors
- Epoxy powder coated frame
- 'C' M-9 size oxygen bottle carrier
- Lightweight composite aluminium/plastic panels, reduces cart noise while in motion and panels can be replaced economically if damaged



## Warranty Information

- 5 Years Frame

Code	Display Name	Total Length	Total Width	Total Height
WAT901180	Emergency Cart (Red) - 5 Drawer + Oxy Carrier + IV Pole + Cardio Board	520 mm	760 mm	1075 mm