

## ROHO - Hybrid Elite Dual Compartment



permobil

**ROHO**

This innovative cushion combines the skin and soft tissue protection and Shape Fitting Technology® of a ROHO DRY FLOATATION® cushion with stability from a contoured foam base. The shape-fitting air insert provides critical pressure redistribution around the pelvis area for those with an obliquity. And the two independent compartments let you set inflation levels to your exact shape and ideal position.

- Combines skin/soft tissue protection and Shape Fitting Technology with stability from a contoured foam base
- Adjustable to shape and weight in the pelvic loading area
- Adjustable to changes in individual's position over time in the pelvic loading area
- Hydrostatic distribution of load in the pelvic loading area
- Provides up to approximately 7 cm of immersion
- Adds little weight to entire mobility system
- Air-filled cushion insert can be easily cleaned and disinfected
- Cover is machine washable and can be disinfected, foam base may be wiped with damp cloth
- Neoprene cushion insert can be patched or repaired to extend use
- Contoured foam base and cushion insert are fluid resistant
- Foam insert is contoured to provide support to the legs
- Foam insert is easily wiped clean and disinfected

### Warranty Information

- Cushion: 2 Years
- Cover: 6 Months

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
PTC594952	ROHO Hybrid Elite - Dual Compartment - 16 x 16in
PTC594953	ROHO Hybrid Elite - Dual Compartment - 16 x 18in
PTC594954	ROHO Hybrid Elite - Dual Compartment - 18 x 18in
PTC594955	ROHO Hybrid Elite - Dual Compartment - 20 x 18in
PTC594956	ROHO Hybrid Elite - Dual Compartment - 15 x 15in
PTC594957	ROHO Hybrid Elite - Dual Compartment - 17 x 16in
PTC594958	ROHO Hybrid Elite - Dual Compartment - 18 x 16in
PTC594959	ROHO Hybrid Elite - Dual Compartment - 18 x 17in
PTC594961	ROHO Hybrid Elite - Dual Compartment - 20 x 20in
PTC594962	ROHO Hybrid Elite - Dual Compartment - 14 x 14in
PTC594963	ROHO Hybrid Elite - Dual Compartment - 16 x 20in
PTC594964	ROHO Hybrid Elite - Dual Compartment - 14 x 15in
PTC594965	ROHO Hybrid Elite - Dual Compartment - 14 x 16in
PTC594966	ROHO Hybrid Elite - Dual Compartment - 15 x 16in
PTC594967	ROHO Hybrid Elite - Dual Compartment - 15 x 17in
PTC594968	ROHO Hybrid Elite - Dual Compartment - 16 x 15in
PTC594969	ROHO Hybrid Elite - Dual Compartment - 16 x 17in

Code	Display Name
PTC594971	ROHO Hybrid Elite - Dual Compartment - 17 x 18in
PTC594972	ROHO Hybrid Elite - Dual Compartment - 18 x 20in
PTC594973	ROHO Hybrid Elite - Dual Compartment - 22 x 18in
PTC594974	ROHO Hybrid Elite - Dual Compartment - 22 x 20in
PTC594975	ROHO Hybrid Elite - Dual Compartment - 24 x 18in
PTC594976	ROHO Hybrid Elite - Dual Compartment - 24 x 20in