

ROHO - Hybrid Elite Dual Compartment




permobil

ROHO

This innovative cushion combines the skin and soft tissue protection and Shape Fitting Technology® of a ROHO DRY FLOATATION® cushion with stability from a contoured foam base. The shape-fitting air insert provides critical pressure redistribution around the pelvis area for those with an obliquity. And the two independent compartments let you set inflation levels to your exact shape and ideal position.

- Combines skin/soft tissue protection and Shape Fitting Technology with stability from a contoured foam base
- Adjustable to shape and weight in the pelvic loading area
- Adjustable to changes in individual's position over time in the pelvic loading area
- Hydrostatic distribution of load in the pelvic loading area
- Provides up to approximately 7 cm of immersion
- Adds little weight to entire mobility system
- Air-filled cushion insert can be easily cleaned and disinfected
- Cover is machine washable and can be disinfected, foam base may be wiped with damp cloth
- Neoprene cushion insert can be patched or repaired to extend use
- Contoured foam base and cushion insert are fluid resistant
- Foam insert is contoured to provide support to the legs
- Foam insert is easily wiped clean and disinfected

 Contracted DVA product Contract approved

Warranty Information

- Cushion: 2 Years
- Cover: 6 Months

| Code | Display Name |
|-----------|--|
| PTC594952 | ROHO Hybrid Elite - Dual Compartment - 16 x 16in |
| PTC594953 | ROHO Hybrid Elite - Dual Compartment - 16 x 18in |
| PTC594954 | ROHO Hybrid Elite - Dual Compartment - 18 x 18in |
| PTC594955 | ROHO Hybrid Elite - Dual Compartment - 20 x 18in |
| PTC594956 | ROHO Hybrid Elite - Dual Compartment - 15 x 15in |
| PTC594957 | ROHO Hybrid Elite - Dual Compartment - 17 x 16in |
| PTC594958 | ROHO Hybrid Elite - Dual Compartment - 18 x 16in |
| PTC594959 | ROHO Hybrid Elite - Dual Compartment - 18 x 17in |
| PTC594961 | ROHO Hybrid Elite - Dual Compartment - 20 x 20in |
| PTC594962 | ROHO Hybrid Elite - Dual Compartment - 14 x 14in |
| PTC594963 | ROHO Hybrid Elite - Dual Compartment - 16 x 20in |
| PTC594964 | ROHO Hybrid Elite - Dual Compartment - 14 x 15in |
| PTC594965 | ROHO Hybrid Elite - Dual Compartment - 14 x 16in |
| PTC594966 | ROHO Hybrid Elite - Dual Compartment - 15 x 16in |
| PTC594967 | ROHO Hybrid Elite - Dual Compartment - 15 x 17in |
| PTC594968 | ROHO Hybrid Elite - Dual Compartment - 16 x 15in |
| PTC594969 | ROHO Hybrid Elite - Dual Compartment - 16 x 17in |

| Code | Display Name |
|-----------|--|
| PTC594971 | ROHO Hybrid Elite - Dual Compartment - 17 x 18in |
| PTC594972 | ROHO Hybrid Elite - Dual Compartment - 18 x 20in |
| PTC594973 | ROHO Hybrid Elite - Dual Compartment - 22 x 18in |
| PTC594974 | ROHO Hybrid Elite - Dual Compartment - 22 x 20in |
| PTC594975 | ROHO Hybrid Elite - Dual Compartment - 24 x 18in |
| PTC594976 | ROHO Hybrid Elite - Dual Compartment - 24 x 20in |