

## Aidacare FL250 Floorline Bed





- Under Bed Clearance While the FL250 achieves a low height of 110 mm, full clearance under the bed is still maintained, minimising the risk of foot or cord impingement
- SmartDescent<sup>™</sup> Automatically slows descent speed over the final 300 mm of lowering to the floor level, allowing possible impingements to be identified and moved
- Hidden Bed Mechanisms: The bed mechanisms are attached under the platform and are not visible. This helps the room feel less clinical
- Low Height 11 cm: Adding a 15 cm mattress gives you a total height below 30 cm thus significantly reducing out-of-bed fall injury risks
- Auto-Regression on back & knee rests: Reduces skin shear forces during operation thus improving clinical outcomes and reducing the need for expensive mattress solutions.
- High SWL 250 kg & base-extension to 100 & 200 mm: The same bed platform can meet a wider profile of residents
- BrakeFirst™ Automatic Safety Brake Activation at lowest height: Reduces the risk of injury during bed exit in the morning as the bed will lock when at its lowest height
- Single Button Seat Position: Easy transfer of weight from lower back to bottom for pressure-care and seated positioning also helps improves lung function for resident
- ExitSafe™ Single Button Ingress/Exit height: Your staff don't have to think about what height to set the bed. Via one button the vast majority of needs are met
- Attendant Control Panel Allows the caregiver to lock out certain functions and also ensures bed operation is not just restricted to the handset
- Grey Colour Scheme & Narrow-profile frame: Keeps the focus on bedding and furnishings/finishings rather than the bed-frame

Contracted DVA product

Ocontract approved

#### **Support @ Home**

Function: Beds and bed equipment

Category: Prescribed

Click **here** for more information

### **Compliance**

IEC 60601-1:2005 IEC 60601-2-52:2005

#### **Accessories**

Click on the cart icon to display the range of optional accessories



- Underbed light standard: Helps to reduce incidence of resident and caregiver falls
- Rubber & Nylon castors: Castors remain on the ground during braking, thus not damaging floors
- King Single Extension Kit (optional retrofit): Allows you to change the profile of your beds (with minimal cost) to meet resident profile changes
- KTL coated anti-corrosion steel: Improved durability for lower lifecycle cost
- Low-maintenance Braking System: Fewer headaches for your MO's!

# Warranty Information

- 10 Years Frame
- 3 Years Electrics
- 1 Year Handset

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
	Aidacare FL250 Floorline Bed - Single (Bed ends sep) - Incl: UBL & Extension	2000 mm	900 mm	110 – 740 mm	250 kg	Sublime Teak
BEB046150	Aidacare FL250 Floorline Bed - King Single (Bed ends sep) - Incl: UBL & Extension	2000 mm	1070 mm	110 - 740 mm	250 kg	Milano Walnut