

## Spex Vigour - Hi Back Support



- Solid full height back support available in three heights for maximum trunk and shoulder support.
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- In combination with lateral supports offers exceptional trunk and thoracic stability
- Approx. Vigour Hi weight is from 1 to 3.3kg

✓ Contract approved

### Support @ Home

**Function:** Accessories for wheelchair seating

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 2 Years

Code	Display Name	Safe Working Load	Backrest Height	Backrest Width
PTB535018	Vigour Kidz - Hi Back Support For 10in W/C (14in High) Rectangular/Full Height Support	90 kg	14 in	9 in
PTB535020	Vigour Kidz - Hi Extended Height Back Support For 10in	90 kg	16 in	9 in
PTB535026	Vigour Kidz - Hi Back Support For 11in W/C (14in High) Rectangular/Full Height Support	90 kg	14 in	10 in
PTB535028	Vigour Kidz - Hi Extended Height Back Support For 11in	90 kg	16 in	10 in
PTB535034	Vigour Kidz - Hi Back Support For 12in W/C (16in High) Rectangular/Full Height Support	90 kg	16 in	11 in
PTB535036	Vigour Kidz - Hi Extended Height Back Support For 12in	90 kg	18 in	11 in
PTB535042	Vigour Kidz - Hi Back Support For 13in W/C (16in High) Rectangular/Full Height Support	90 kg	16 in	12 in

Code	Display Name	Safe Working Load	Backrest Height	Backrest Width
PTB535044	Vigour Kidz - Hi Extended Height Back Support For 13in	90 kg	18 in	12 in
PTB535050	Vigour - Hi Back Support For 14in W/C (18in High) Rectangular/Full Height Support	90 kg	18 in	13 in
PTB535052	Vigour - Hi Extended Height Back Support For 14in W/C (20in High) Rectangular/Extra Height Support	90 kg	20 in	13 in
PTB535058	Vigour - Hi Back Support For 15in W/C (18in High) Rectangular/Full Height Support	90 kg	18 in	14 in
PTB535060	Vigour - Hi Extended Height Back Support For 15in W/C (20in High) Rectangular/Extra Height Support	90 kg	20 in	14 in
PTB535066	Vigour - Hi Back Support For 16in W/C (20in High) Rectangular/Full Height Support	160 kg	20 in	15 in
PTB535068	Vigour - Hi Extended Height Back Support For 16in W/C (23in High) Rectangular/Extra Height Support	160 kg	23 in	15 in
PTB535074	Vigour - Hi Back Support For 17in W/C (20in High) Rectangular/Full Height Support	160 kg	20 in	16 in
PTB535076	Vigour - Hi Extended Height Back Support For 17in W/C (23in High) Rectangular/Extra Height Support	160 kg	23 in	16 in
PTB535082	Vigour - Hi Back Support For 18in W/C (20in High) Rectangular/Full Height Support	160 kg	20 in	17 in
PTB535084	Vigour - Hi Extended Height Back Support For 18in W/C (23in High) Rectangular/Extra Height Support	160 kg	23 in	17 in
PTB535090	Vigour - Hi Back Support For 19in W/C (20in High) Rectangular/Full Height Support	160 kg	20 in	18 in
PTB535092	Vigour - Hi Extended Height Back Support For 19in W/C (23in High) Rectangular/Extra Height Support	160 kg	23 in	18 in
PTB535098	Vigour - Hi Back Support For 20in W/C (20in High) Rectangular/Full Height Support	160 kg	20 in	19 in
PTB535100	Vigour - Hi Extended Height Back Support For 20in W/C (23in High) Rectangular/Extra Height Support	160 kg	23 in	19 in