

## FreeStyle Libre 2 Continuous Blood Glucose Monitor - Sensor



The FreeStyle Libre 2 sensor is factory calibrated and the results are stable for the life of the sensor— so you don't need to calibrate it again.

- It continuously measures the glucose concentration in your interstitial fluid, and stores 8 hours of data.
- Optional alarms - With optional real-time glucose alarms, rest assured you'll know when your glucose is too high or too low.
- Alarm - Receive an alarm when your glucose is too low or too high.
- Check - Check the FreeStyle LibreLink app, or your FreeStyle Libre 2 reader to see your glucose reading, trend arrow, and 8-hour history.
- Act - Use the information to take appropriate action as needed.
- Small and discreet - it is designed to be placed on the back of your upper arm and can stay on for up to 14 days.
- Easy to use<sup>3</sup> & comfortable to wear - painless to apply and easy to wear.
- Swim, shower or exercise — the sensor is designed to be water-resistant and worn while bathing, showering, swimming and exercising.

✓ Contracted DVA product

