


Aspire ComfiMotion Fixed Allure Headboard



Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs—improving comfort for users and making life easier for carers.

Fixed headboards attach directly to the bed frame thereby moving with the bed during tilt and height adjustments.

- By moving with the bed this reduces pillow displacement and supports the pillow and head position during tilt
- Pivoting design enhances wall protection by pivoting upon wall contact
- The Bed/Headboard can be positioned in direct proximity to a wall and still perform head-down tilt function
- Built-in handset loops (Double/Queen) keep controls within convenient reach
- Widths to match bed sizes, all 680mm high and 120mm deep

 Contract approved

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Colour
BEA830154	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - Long Double	749	1	120 mm	1370 mm	680 mm	-
BEA830156	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - Queen	799	1	120 mm	1520 mm	680 mm	Portobello Grey
BEA830152	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - King Single	649	1	120 mm	1070 mm	680 mm	Portobello Grey

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Colour
BEA830150	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - Long Single	599	1	120 mm	920 mm	680 mm	Portobello Grey