

Mangar Elk - Lift Cushion



The ELK is a lifting cushion designed to lift from the floor. The small footprint design allows for the ELK to be used in confined spaces, where traditional lifting equipment may not fit. Ideal for use within a home environment.

- Inflate effortlessly using the Airflo Compressor with just a press of a button. No complex setups or extensive training is needed – it is that simple! Whether indoors or outdoors, this system is designed for seamless operation by a single caregiver.
- Easy to clean using a proprietary, non-abrasive mild liquid cleaner, simply deflate, fold, and stow away in the provided carry bag – ideal for infection control management.
- Capable of lifting individuals weighing up to 450kg from the floor, it provides unparalleled stability and rigidity. Ideal for carers tending to individuals with neurological disorders prone to falls, it offers peace of mind in every lift.
- Experience a lift in both safety and dignity. By reducing the risk of musculoskeletal injuries for caregivers, it is a game-changer in care routines.



Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name
LSS403820	Mangar Elk - Lift Cushion - 450kg