

Megan Lowback Chair



- Contemporary and stylish accent chair suited to modern interiors
- A unique take on a classic tub chair with flared arms and a low back
- Versatile piece that can be used in a bedroom, lounge or waiting area
- Mid-century inspired legs with a warm teak finish
- Firm, yet soft cushioning provides ample support
- Fabric colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL052160	Megan Lowback Chair	710mm	550mm	735mm