

## Aspire ComfiMotion Care Bed Package - Chesterfield



The Aspire ComfiMotion Care Bed Package brings together advanced comfort, safety and style in a complete solution for home or care environments. Supplied with the Chesterfield Free Standing Headboard as standard, this package combines the functionality of the ComfiMotion Care Bed with a modern, stylish finish - delivering both clinical support and a welcoming bedroom aesthetic.



- Night Assist: Luminous underbed halo for safe entry and exit throughout the night
- Press and hold safety handset: Intuitive, large button interface enhances usability
- Enclosed frame design: Ensures optimal positioning when getting into the adjustable bed and safely covers the moving base platform
- Safe-Lowering Technology: Ensures that back, leg and Hi-Lo functions only lower under the force of gravity reducing potential for accidents or injuries
- Box-on-box design facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- German engineered motors for smooth, quiet adjustment
- \*Total Height with:
  - Low Profile Fixed Legs 280 – 580 mm (available upon request)
  - 75mm Standard Castors 350 - 650 mm
  - 125mm Accessibility Castors 400 - 700 mm (optional)
- Mattress choice may increase or decrease overall bed height
- Mattress not included
- **Prescribing Therapists:** Click here to download our Clinical Assistance Guide

### Warranty Information

- 10 Years Frame
- 2 Years Electronics

- 1 Year Handset, Cables & Fabric
- 1 Year Headboard

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour
BEB046720V2	Aspire ComfiMotion Care Bed Package - Chesterfield - Long Single (inc. Freestanding Headboard)	2030 mm	920 mm	280 - 580 mm	180 kg	90 kg	Portobello Grey
BEB046730V2	Aspire ComfiMotion Care Bed Package - Chesterfield - King Single (inc. Freestanding Headboard)	2030 mm	1070 mm	280 - 580 mm	180 kg	100 kg	Portobello Grey