

## Paediatric Mobiliser



- Mobilisor range is a sturdy, folding, height - adjustable walking frame for children
- Includes fixed front castors and rear glides
- Width: 440 mm (Child), 495 mm (Junior)
- Handle Height: 600 - 790 mm (Child), 670 - 960 mm (Junior)



### Compliance

ISO 11199-2:2005

## Warranty Information

- 1 Year

| Code      | Display Name       | Safe Working Load |
|-----------|--------------------|-------------------|
| WAF704200 | Mobiliser - Child  | 60 kg             |
| WAF708600 | Mobiliser - Junior | 80 kg             |