

Aspire Homecare Shower Stool



- Allows for a comfortable bathing experience for users who require that little bit of extra support
- Durable high grip feet for stability and safety
- Adjustable legs to suit a wide range of users
- Raised armrests assist the user in transferring on and off
- Width between the arms: 465 mm

- ✔ Contracted DVA product
- ✔ Contract approved

Warranty Information

- 1 Year

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information



Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS119035	Aspire Homecare Shower Stool	525 mm	590 - 720 mm	100 kg	340 mm	410 mm	440 - 570 mm	440 - 570 mm