

TabTimer ADVANCE Careousel Automatic Pill Dispenser



TabTimer Careousel Automatic Pill Dispenser



Advanced Technology, Enhanced Security - The Advance medication dispenser incorporates a number of important and advanced features designed to help both the user remain compliant with their medication regime and further support their family members, carers, and pharmacists. The alarm can be programmed to ring for up to 5 hours. With 28 sections and a possible 24 alarms per day, the Pivotell Advance automatic pill dispenser can accommodate the most complicated of medication regimes.

- **1-28 Daily Alarm Times** - Accommodates the most complicated of medication regimes
- **29 Compartments** - For easy pharmacy filling (one compartment is left empty)
- **Dual / Multi Dose programming** - Two or more compartments can be programmed to alarm at the same time. For those with many or large pills
- **Start delay function** - Delays activation to allow pharmacists or family members to easily fill in advance away from the user's home
- **Pharmacy label carrier and label case** - supplied to comply with pharmacy labelling requirements
- **Upside down feature** - The medication tray rotation is disabled if the dispenser is left upside down, to avoid pill spillage. The alarm will still sound at a dosage time.
- **Day / time discs supplied** for 1 / 2 / 3 / 4 times daily dosage, and one blank disc for five or more doses per day
- **Battery level meter** and warning indicator
- **Four different alarm sounds**
- **Flashing light** for those with hearing impairment
- **Security** - a self-locking shutter over the lid opening prevents access to the

medication tray except at the alarm time. This prevents double dosing or forced movement of the tray by the user.

- **Display** - LCD shows date, time, next dose, doses remaining, and battery level.
- Same diameter and slightly deeper than Pivotell® Mk3/11 Dispenser
- Incorporates antibacterial silver ion technology for microbial protection

PHARMACY REFILLING? - If the patient will be using a pharmacy to refill the device it is highly advisable that you purchase the bundle with spare tray and tamper-evident stickers to assist the pharmacy with ease of refilling. The tamper-evident stickers help pharmacies comply with the QCPP Quality Care Pharmacy Program requirements. The TabTimer Automatic Pill Dispenser meets the requirements of the QCPP (Quality Care Pharmacy Program) for a dose administration aid, when used with the supplied medicine card and optional tamper-evident stickers.

Code	Display Name
DLG275010	TabTimer ADVANCE Careousel Automatic Pill Dispenser - Solid Lid