

## BODY ASSIST 410S THERMAL KNEE SUPPORT



- Therapeutic warmth, superior compression and protection.
- Contoured fit - closed patella.
- CR neoprene is 100% latex free lightweight medical grade neoprene.
- Fits left or right knee.
- Washable, hygienic, reversible and lined both sides.
- Heavy duty outer lining for durability.
- Medoxylin liner provides thermal insulation, moisture draw capacity and air circulation to keep skin well oxygenated and comfortable. Application:
- For minor support helping to relieve pain and promote healing of tendonitis, bursitis, arthritis and other nonspecific knee conditions.

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Lower limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS307720	BODY ASSIST 410S THERMAL KNEE SUPPORT - SML
FBS307730	BODY ASSIST 410S THERMAL KNEE SUPPORT - MED
FBS307740	BODY ASSIST 410S THERMAL KNEE SUPPORT - LGE
FBS307750	BODY ASSIST 410S THERMAL KNEE SUPPORT - XL
FBS307760	BODY ASSIST 410S THERMAL KNEE SUPPORT - XXL