

Aspire Hybrid Patient Sling



- Hybrid sling facilitates both hygiene and general purpose manual handling
- Breathable mesh promotes air circulation, reducing moisture and heat build-up while helping to maintain skin integrity
- Intuitive shape allows ease of access for toileting and dressing
- Security belt and padded leg straps support the trunk and thighs
- Full coverage shoulder and head support
- Durable loop webbing designed not to tangle or ‘cork-screw’ during cleaning
- Embroidered identification labelling ensures equipment traceability, supporting inventory control and maintenance protocols
- Machine washable up to 85°C



Warranty Information

- 1 Year



Code	Display Name	Safe Working Load
LSS390856	Aspire Hybrid Patient Sling Loops - X Large	300 kg
LSS390854	Aspire Hybrid Patient Sling Loops - Large	300 kg
LSS390852	Aspire Hybrid Patient Sling Loops - Medium	200 kg
LSS390850	Aspire Hybrid Patient Sling Loops - Small	200 kg