

## MoliCare Premium Mobile



- MoliCare Premium Mobile pants look and feel like regular underwear through anatomical shape, soft material and Lycra for a comfortable fit and discreetness
- Curly fibre layer provides a highly antibacterial and pH 5.5 neutral effect to maintain healthy skin
- MoliCare Premium Mobile pants are dermatologically tested
- The absorbent core neutralises odour and holds liquids quickly for a dry feeling
- Maximum leakage protection is secured based on high-quality material and the two-part cuff system
- Tear-open side seams assure a convenient removal
- Coloured Lycra indicate size, absorbency and front garment



### Compliance

No applicable compliance standard

Code	Display Name	Capacity
CON253320	MoliCare Mobile 6D - Small (Pkt 14)	6 Drops (1361 - 2140ml)
CON253325	MoliCare Mobile 6D - Medium (Pkt 14)	6 Drops (1361 - 2140ml)
CON253330	MoliCare Mobile 6D - XL (Pkt 14)	6 Drops (1361 - 2140ml)
CON253335	MoliCare Mobile 8D - Small (Pkt 14)	8 Drops (1791 - 2421ml)
CON253341	MoliCare Mobile 8D - Medium (Pkt 14)	8 Drops (1791 - 2421ml)
CON253345	MoliCare Mobile 8D - Large (Pkt 14)	8 Drops (1791 - 2421ml)
CON254600	MoliCare Mobile 6D - Large (Pkt 14)	6 Drops (1361 - 2140ml)

Code	Display Name	Capacity
CON254700	MoliCare Mobile 8D - XS (Pkt 14)	8 Drops (1791 - 2421ml)
CON254710	MoliCare Mobile 8D - XL (Pkt 14)	8 Drops (1791 - 2421ml)
CON254715	MOLICARE Mobile 10D - Medium (Pkt 14)	10 Drops (2626 - 2757ml)
CON254720	MOLICARE Mobile 10D - Large (Pkt 14)	10 Drops (2626 - 2757ml)