

Tray Clearing Trolley



- Tray Clearing Trolley is ideal for hospital or nursing home environment
- Stainless steel construction featuring two or three shelves design
- Four 150 mm swivel castors

Warranty Information

- 3 Years



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height	Material
WAT824300	Tray Clearing Trolley 2 Shelf with Receptacles	500 mm	1100 mm	1000 mm	Stainless Steel
WAT824500	Tray Clearing Trolley 3 Shelf with Receptacles	500 mm	810 mm	965 mm	Stainless Steel