

Tray Clearing Trolley



- Tray Clearing Trolley is ideal for hospital or nursing home environment
- Stainless steel construction featuring two or three shelves design
- Four 150 mm swivel castors



Compliance

No applicable compliance standard

Warranty Information

- 3 Years

Code	Display Name	Total Length	Total Width	Total Height
WAT824300	Tray Clearing Trolley 2 Shelf with Receptacles	500 mm	1100 mm	1000 mm
WAT824500	Tray Clearing Trolley 3 Shelf with Receptacles	500 mm	810 mm	965 mm