

Mini Duke Bedside Cabinet



- Classic bedside cabinet to suit contemporary interiors
- Three drawers offer plenty of storage for personal items
- Ample space for a bedside lamp, glasses and book
- Lockable top drawer keeps valuables safe
- Raised legs keeps the base off the ground and enable cleaning underneath
- Smooth slimline handles are easy to hold
- Dimensions: 700 mm (H) x 500 mm (L) x 480 mm (W)
- Pictured in Milano Walnut, colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height	Colour
BEA826200	Mini Duke Bedside Cabinet - Laminate	450 mm	500 mm	700 mm	-
BEA826210	Mini Duke Bedside Cabinet - Timber	450 mm	500 mm	700 mm	-
BEA826201	Duke Mini Bedside Cabinet - Laminate/Melamine - Sublime Teak	450 mm	500 mm	700 mm	Subline Teak
BEA826202	Duke Mini Bedside Cabinet - Laminate/Melamine - Milano Walnut	450 mm	500 mm	700 mm	Milano Walnut