

Body Assist 43H Hinged Knee Brace



- Rigid hinges minimise excessive medio - lateral movement of the knee joint, while offering excellent support & stability
- Made from lightweight medical grade orthopaedic elastic & features front closures for easier application
- Top & bottom straps lock the brace in position & minimise movement

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS280770	BodyAssist - Knee Hinged - Elastic Wrap - O/Patella - Sml	up to 320 mm
FBS284300	BodyAssist - Knee Hinged - Elastic Wrap - O/Patella - Med	320 - 370 mm
FBS284310	BodyAssist - Knee Hinged - Elastic Wrap - O/Patella - Lge	370 - 420 mm
FBS284315	BodyAssist - Knee Hinged - Elastic Wrap - O/Patella - XL	420 - 470 mm
FBS284220	BodyAssist - Knee Hinged - Elastic Wrap - O/Patella - XXL	470 mm+