

Kerry General Purpose Sling with Head Support - Polyester



- Attaches to hoist to enable a carer to lift and manoeuvre a less mobile user
- Provides excellent back and head support for patients requiring support during transfers
- Quick and easy to fit
- Ideal for non weight bearing patients
- Colour coded attachment and positioning straps
- Machine washable
- FIM Scores: 1 and 2



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

| Code | Display Name | Safe Working Load | Material |
|-----------|---|-------------------|-----------|
| LSS381500 | General Purpose Sling with Head Support - Small | 250 kg | Polyester |
| LSS381700 | General Purpose Sling with Head Support - Large | 250 kg | Polyester |
| LSS381750 | General Purpose Sling with Head Support - Extra Large | 250 kg | Polyester |