

Aspire Overbed Table Frame




- Stable and secure table top to conveniently accommodate meals, drinks, reading materials, tablet devices and other items within easy arms reach from a bed or chair
- Height adjustment lever is easily accessible to users and carers
- Height fixing knob will lock the height and prevent adjustment by users
- Two locking castors for securing table position by the carer to prevent it rolling away

Warranty Information

- 2 Years



 Contract approved

ndis

Item Code:

05_180315322_0103_1_2

Minor Bed – Composite:

05_428_0103_1_2

Click [here](#) for more information

Compliance

ISO7172

| Code | Display Name | RRP | Box Quantity | Total Height |
|-----------|--|-----|--------------|---------------|
| BEA006685 | Aspire Overbed Table - Laminated & Customs - Frame - Gunmetal Grey | 236 | 1 | 730 - 1095 mm |
| BEA006690 | Aspire Overbed Table - Laminated & Customs - Frame - Cream | 236 | 1 | 730 - 1095 mm |