

Queens Angled Cutlery



- A range of stainless steel utensils with an angled head mounted in the Queens built-up handle, enabling those with limited wrist movement to eat more easily
- The knife has a blade bent at an angle to the handle and can be used with a rocker action
- The fork and spoon are available as right or left handed versions
- Handle Length: 98 mm
- Handle Diameter: 30 mm
- Product Weight: 45 - 55 g

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Eating and drinking

Category: Low risk

Click [here](#) for more information

NDIS

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Unit Weight
DLK275820	Angled Knife	45 g
DLK275830	Angled Fork - Left	55 g
DLK275840	Angled Fork - Right	55 g

Code	Display Name	Unit Weight
DLK275850	Angled Spoon - Left	55 g
DLK275860	Angled Spoon - Right	55 g