

Aspire Clamp-On Bed Support Rail



- Designed with clinical functionality in mind, the Aspire Clamp-On Bed Support Rail enhances patient safety during repositioning, transfers, and transitions to and from mobility aids like wheelchairs
- The clamp design ensures consistent setup throughout different facilities, with simple and seamless training and deployment
- Designed to fit the Aspire Lifestyle Community and AC4 Beds



Warranty Information

- 5 Years

Code	Display Name	RRP	Box Quantity	Total Width	Total Height
BEA829105	Aspire Clamp-On Bed Support Rail	129	1	385 mm	425 mm