

Overbed Table - Split Table Top



- The split top can be tilted and fixed at 5 different positions to suit the user while smaller surface remains horizontal to hold drinks or other items
- Locking tilt mechanism eliminates surface movement in all positions
- Wide base model fits around most lift recliners and chairs
- Tilting top can be positioned on either left side or right side of fixed top as desired
- Additional Specs:
 - Table Top Width: 430 mm
 - Table Top Length: 800 mm
 - Tilt Section Length: 600 mm
 - Fixed Section Width: 200 mm
 - Internal Base Width: 800 mm

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Tables

Category: Prescribed

Click [here](#) for more information



Warranty Information

- 1 Year

Code	Display Name	Total Height
BEA028110	Overbed Table - Split Table Top	750 - 1100 mm