

Jobst Bella Lite - Standard Arm Sleeve (20 - 30 mmHg)



- Soft and air-permeable
- Lightweight knit construction
- Fits a wide range of hand and arm shapes
- Fits both the left or right arm/hand
- The combination sleeve/glove can be worn separately together

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information



Code	Display Name	Calf/Arm Circumference	Elbow Circumference	Wrist Circumference
FBS303478	Jobst Bella Lite - Ready To Wear - Standard 20-30mmHg Sml Long	220 - 400 mm	180 - 230 mm	140 - 200 mm
FBS303479	Jobst Bella Lite - Ready To Wear - Standard 20-30mmHg Sml Reg	220 - 400 mm	180 - 230 mm	140 - 200 mm
FBS303481	Jobst Bella Lite - Ready To Wear - Standard 20-30mmHg Med Reg	270 - 450 mm	230 - 290 mm	160 - 220 mm
FBS303483	Jobst Bella Lite - Ready To Wear - Standard 20-30mmHg Lge Reg	330 - 510 mm	290 - 340 mm	180 - 240 mm