

DONJOY 82-9613 DUOFORM WRIST SUPPORT



- Easy application
- Universal model, fits either right or left hand
- Repositionable dorsal stay, provides appropriate tailored support
- Ergonomic palmar splint to control wrist flexion
- Removable modular components for universal use
- Wide opening for easy application
- Three straps for adjustable compression
- Breathable 3D material
- Comfortable, flexible, lightweight design Application:
- Wrist sprains and trauma
- Pre and post-operative immobilization

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS305800	DONJOY 82-9613 DUOFORM WRIST SUPPORT - XS
FBS305810	DONJOY 82-9613 DUOFORM WRIST SUPPORT - SML
FBS305820	DONJOY 82-9613 DUOFORM WRIST SUPPORT - MED
FBS305830	DONJOY 82-9613 DUOFORM WRIST SUPPORT - LGE
FBS305840	DONJOY 82-9613 DUOFORM WRIST SUPPORT - XL