

Donjoy Ultra Sling III



- Moisture-wicking mesh fabric for greater comfort, particularly in longer-term use
- Helps prevent post-operative internal rotation contractures
- Promotes auxiliary air exchange to reduce risk of secondary infections and prevents internal rotation contractures
- Easy-open front panel encourages forearm exercises
- Allows arm to be positioned in a variety of positions dependent on post-operative requirements (either more posterior or anterior)
- The De-Rotation Strap hooks to the UltraSling III and prevents internal rotation by securely holding the patient's arm in a neutral position
- Quick-release buckles make it easy to take on and off for exercise/therapy

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Wrist Circumference
FBS301370	Donjoy Ultra Sling III - Medium	140 - 230 mm
FBS301371	Sling - Donjoy - Ultra Sling III - Large	140 - 230 mm