

## Aspire ComfiMotion Freestanding Chesterfield Headboard



Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs - improving comfort for users and making life easier for carers.

Freestanding headboards stand independently behind the bed and can be moved or removed as needed. Contemporary in design, the headboard remains stationary while the bed is raised or lowered, not impacting any wall decoration above the headboard.

- Not attached to the bed - offers full positioning flexibility
- Quick to install or remove without tools
- Softened surface that supports pillows and protects the head
- Enhances room style without limiting bed function
- Supports high-needs equipment integration in care environments
- Available in two modern styles, both at 1140mm high and 110mm deep

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

IEC 60601-2-52:2009

| Code      | Display Name   |
|-----------|--|
| BEA830100 | Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - Long Single |
| BEA830110 | Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - King Single |
| BEA830112 | Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - Long Double |
| BEA830114 | Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - Queen       |
| BEA830116 | Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - King        |

