

Oxford Ascend Standing Lifter



- With it's great range of movement, the Ascend can transfer a patient from a low-seated position to fully extended with ease
- The Ascend is compact and easy to manoeuvre in confined spaces such as bathrooms and has many key design features that make it a truly flexible piece of equipment
- Its active lifting motion has been designed to encourage user participation, which ultimately promotes patient independence and well-being
- The kneepad can be adjusted within a 150 mm operating range to improve patient comfort during repositioning
- The removable foot tray allows the Ascend to be used as a walking aid - great for patients that require gait-training
- Ergonomic foot pad to assist with initiating movement
- Careful attention to product design ensures correct ergonomic usage, protecting both the patient and carer from injury at all times
- Item Weight: 47 kg

✔ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Warranty Information

- 2 Years

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
LSS377350	Ascend Patient Lifter	1100 mm	630 mm	1070 - 1710 mm	170 kg	Aluminium	47 kg	990 - 1695 mm	510 - 885 mm	True