

## Oxford Ascend Standing Lifter



- With it's great range of movement, the Ascend can transfer a patient from a low-seated position to fully extended with ease
- The Ascend is compact and easy to manoeuvre in confined spaces such as bathrooms and has many key design features that make it a truly flexible piece of equipment
- Its active lifting motion has been designed to encourage user participation, which ultimately promotes patient independence and well-being
- The kneepad can be adjusted within a 150 mm operating range to improve patient comfort during repositioning
- The removable foot tray allows the Ascend to be used as a walking aid - great for patients that require gait-training
- Ergonomic foot pad to assist with initiating movement
- Careful attention to product design ensures correct ergonomic usage, protecting both the patient and carer from injury at all times
- Item Weight: 47 kg



### ndis

#### Item Code:

05\_123603161\_0105\_1\_2

#### Hoist –

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

## Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS377350	Ascend Patient Lifter	1100 mm	630 mm	1070 - 1710 mm	990 - 1695 mm	510 - 885 mm	True	170 kg	Aluminium