

Aspire ComfiMotion Activ Care Bed





The Activ Care bed combines the best of in-home care functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

- Ultra Low Height: Minimum platform height of 195 mm reduces falls risk, improving user safety
- Customisable Safe Entry/Exit Height: Programmable single button press ensures perfect bed height to maximise independence
- Under Bed Clearance: provides clear access for patient lifters and overbed tables
- Battery Backup: Provides confidence and maintains independence during power outages
- Night Assist: Luminous underbed halo for safe entry and exit
- Illuminated Handset: Large button interface is intuitive and highly visible at all times
- Massage: Relaxing vibration modes for personalised neck, shoulder and back relief
- Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks
- Head Tilt: Relaxes neck muscles, providing optimal support whilst reading or watching TV
- Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea
- Auto–Regression in the Backrest: Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution
- **Clinical Positioning:** One touch raises back and legs to alleviate cardiac exertion, improving resting comfort
- Leg Raise: Promotes circulation, relieving tired, swollen/ fatigued or painful



Compliance

AS/NZS 60601-2-52:2009



limbs

- **Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button
- Zero Gravity: Helps to relax the circulatory system and muscles of the lower back
- Seated Position: Assists with bed entry and exit and upright dining
- Trendelenburg (and Reverse): Can assist care givers when repositioning a user*
- Lockable Castors: Eight heavy duty dual castors
- Headboard and accessories sold separately

Warranty Information

10 Years Limited

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
BEB046745	Aspire ComfiMotion Activ Care Bed - Portobello Grey - Long Single (Headboard Sep)	2030mm	920mm	200 - 775mm 2	50kg	Portobello Grey
BEB046755	Aspire ComfiMotion Activ Care Bed - Portobello Grey - King Single (Headboard Sep)	2030mm	1070mm	200 - 775mm 2	50kg	Portobello Grey
BEB046765	Aspire ComfiMotion Activ Care Bed - Portobello Grey - Long Double (Headboard Sep)	2030mm	1370mm	200 - 675mm 3	60kg	Portobello Grey

^{*}Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
BEB046775	Aspire ComfiMotion Activ Care Bed - Portobello Grey - Queen (Headboard Sep)	2030mm	1520mm	200 - 675mm 3	60kg	Portobello Grey