

## NuboMed X4000 Advanced Mobile Nursing Cart



NuboMed

- Right side facing drawers, duo handles for side operations to maximize storage offerings.
- Hidden and expandable work surface for maximum usage space.
- A variety of accessories and drawer combinations.
- Mechanical lock for security.
- Embedded computer, HD display, design with care.
- Hot-swappable main and auxiliary batteries, 24h uninterrupted usage.
- Imported Panasonic battery core from Japan, safety guaranteed.
- Wireless charging (Phone/PDA), RFID user log-in.
- Ideal For: Daily Nursing/Data Recording



### Warranty Information

- 2 Years

Code	Display Name
WAT908175	NuboMed - X4000 - Advanced Mobile Nursing Cart