

NuboMed E3000 Height Adjustable Cart



NuboMed

- **Material:** Constructed from aluminium alloy
- **Integrated Design:** Made from a single piece of die-cast aluminium alloy, this design integrates the handle and keyboard tray for a screwless, gap-free structure, while also featuring a 28mm raised side guard and an internal battery
- **Height Adjustable:** Pneumatic adjustment by gas spring with a range of 0-300mm
- **LED Light:** Infrared proximity-sensing LED lighting module positioned above keyboard
- **Mouse Pad:** Concealable and rotating (136 x 138mm), can be installed on the left or right
- **Monitor Mount:** Height adjustable from 0-150mm, 15 degrees downward and upward tilt, and portrait-to-landscape rotation
- **Display Panel:** LED status indicators for real-time battery levels and abnormal conditions
- **Button:** PC power on/off button and a screen backlight on/off switch
- **Base:** Designed to prevent tripping and made from a single-piece of aluminium alloy with an antibacterial powder coating finish
- **Castors:** Equipped with 2x 4" lockable and 2x 4" tracking castors
- **Ideal For:** Doctor's Rounds/Training/Conferencing, Daily Nursing



Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Unit Weight
WAT908210	NuboMed E3000 Height Adjustable Cart	653 mm	462 mm	1680 mm	25.4 kg