

Post Operative Steps



- Single and Double steps for bed access
- Single step - 240 mm. Double step - 220 mm and 400 mm
- Black, non-slip ribbed rubber surface
- Rubber stoppers for stability

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Material
BEA014900	Post Op Step - Single	Steel, Rubber
BEA014800	Post Op Step - Double	Steel, Rubber