

## Aspire Male Urinal



Aspire Male Urinals are fully autoclavable and made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Collecting urine and faeces

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 1 Year



Code	Display Name	Capacity
BTT145700	Male Urinal - Autoclavable - Green	1000 ml