

Meiko Internal Rack



- Holder for care utensils and urine bottles
- Stainless steel
- Quick-change

✔ Contract approved

Support @ Home

Function: Collecting urine and faeces

Category: Prescribed

Click [here](#) for more information



Code	Display Name
UPA400400	Meiko Internal Rack - SUH1 - (non-standard in Aust)
UPA400410	Meiko Internal Rack - SUH2 - Quick-change rack (standard in Aust)
UPA400420	Meiko Internal Rack - TopCut B - Quick Change Rack with Bag Empty
UPA400430	Meiko Internal Rack - TopCut H - Bag Emptying Rack