

Kerry Full Body Sling - Mesh



- Attaches to hoist to enable a carer to lift and manoeuvre a less mobile user
- Made from a lightweight mesh which is ideal for showering or bathing and dries quickly
- Divided leg is ideal for toileting patients who require full body support
- Provides excellent back and head support for patients requiring support during transfers
- Colour coded attachment and positioning straps
- Machine washable



ndis

Item Code:

05_123603165_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS382600	Full Body Sling - Small	250 kg	Mesh
LSS382700	Full Body Sling - Medium	250 kg	Mesh
LSS382800	Full Body Sling - Large	250 kg	Mesh