

## Dementia Care Glowing Tape



- A glow in the dark tape containing a lumescence substance which is designed to be visible for more than 100 hours after the lights go out
- It is charged from overhead and natural lighting
- Suitable for adhesion to smooth, hard, non-porous surfaces in dry areas (indoor)
- Self adhesive

### Warranty Information

- 1 Year

- ✔ Contracted DVA product
- ✔ Contract approved

### Support @ Home

**Function:** Orientation

**Category:** Low risk

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width
DCO100125	Glow tape - narrow 25mm W x 10M long roll	10 m	25 mm
DCO100130	Glow tape - wide 50mm W x 10M long roll	10 m	50 mm