

TENA Pad Lady



- Provides an excellent level of security, comfort and discretion
- Best suited to ambulant people or those who are able to toilet themselves
- To be used in conjunction with TENA Fix Stretch Pants or TENA Fix Cotton Special as part of a two-piece system

Warranty Information

- 1 Year



ndis

Item Code:

03_092489060_0103_1_1

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Capacity
CON243400	TENA Pad Lady Extra (Pkt 20)	295 ml
CON243450	TENA Pad Lady Extra Plus (Pkt 8)	350 ml
CON243500	TENA Pad Lady Super (Pkt 30)	515 ml