

Aspire ComfiMotion Bed Collection - Bed Rail



- Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility
- Designed to be used in conjunction with ComfiMotion bed collection - please read variant names for compatibility
- Extendable Version (BEA830205): Ability to extend to prevent falls and contract to allow safe bed egress and ingress (2nd image)
- Additional specifications
 - Height of rail from mattress deck: 565 mm
 - Width of rail: 685 mm

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name
BEA830210	Aspire ComfiMotion CARE Bed - Bed Rail