

## Aspire Mobility Accessory Bags



### Stylish & practical bags for walkers and wheelchairs.



Keep personal items close and secure when you're on the move. Designed for use with most seat walkers and wheelchairs, these bags are ideal for storing everyday essentials like phones, wallets, water bottles or tissues. No need to rely on carers or balance items on your lap.

- **Practical & Spacious:** Plenty of room for daily must-haves like keys, phones, glasses, and snacks
- **Secure & Stable:** Button-press lock allows personal belongings to be stored safely
- **Versatile & Stylish:** Choose from four different designs
- **Universal Fit:** Easily attaches to the front bar of most walkers or rear handles of most wheelchairs

### Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Safe Working Load	Colour
WAA696500	Mobility Accessory Bag - Walkers/Wheelchairs - Leopard	49	1	5 kg	Leopard
WAA696510	Mobility Accessory Bag - Walkers/Wheelchairs - Floral	49	1	5 kg	Floral
WAA696520	Mobility Accessory Bag - Walkers/Wheelchairs - Plaid	49	1	5 kg	Plaid
WAA696530	Mobility Accessory Bag - Walkers/Wheelchairs - Blue	49	1	5 kg	Blue

