

## Kerry Multi Lift Patient Lifter



- Electronic handset control of lifting arm height, single speed with soft start/stop
- Electronic handset control of leg gating
- Weight overload control and indication
- Emergency stop button
- Hand-operated emergency lowering
- Wide, lateral push handle and large diameter rear castors for excellent manoeuvrability
- Small front castors for access under lower furniture
- Rechargeable battery, with low battery warning alarm, that can be removed and charged
- Monitoring and indication of battery charge level
- Maintenance Required indication at set service intervals
- Moulded ABS cover for durability and easy cleaning
- Long-lasting LINAK lifting actuator for lower ongoing maintenance costs



### Warranty Information

- 2 Years