

Caring Mug



- The Caring Mug features a wide stable base with two large contoured and angled handles
- It is designed to be held in a comfortable position for drinking and is useful for people with tremor or who cannot sit upright. The mug is microwaveable
- Weight: 176 g

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Capacity
DLK276510	Clear Mug	300 ml